

FlowZig

THE WIRELESS UROFLOWMETRY



FlowZig

Wireless technology for precise diagnostics

Wireless Uroflowmetry with

- BLUETOOTH® protocol
- Low power consumption
- Band 2400-2485 MHz
- Max. operative distance 5 m
- Compatible with any Laptop or Desktop working with Windows Operating System
- Software included

Analysis




- Uroflowmetry
- Manual and Automatic Analysis
- Siroky maximum and average flow rate nomograms
- Liverpool Nomograms

* Technical data can be changed without any notice



Flowmeter

- Load cell flow trasducer at high resolution and long term accuracy
- Flow range 0 to 100 ml/s
- Accuracy 1% F.S.
- End scale: automatic or manual
- Sampling rate: 10Hz
- A/D converter: 16 bit
- Volume range: 2000 ml



MANUFACTURER

MEDICA S.p.A. Via degli Artigiani, 7 | 41036, Medolla (Modena) Italia

Tel: +39 0535 51159 | Fax: +39 0535 52605 | Email: info@medica-spa.com | www.medica-spa.com

MEDICA

SEAMLESS INNOVATION FOR LIFE QUALITY



CE 0476 UNI EN ISO 9001:2015 UNI CEI EN ISO 13485:2021