

HEMOFILTERS

CRRT FROM NEWBORNS
TO ADULTS



FEMALE LUER LOCK CONNECTION

	DP07HE	DP09HE	DP12HE	DP15HE	DP03HE	D050	D150
Intended use	Adult				Pediatric		
Hollow fiber material	HIGH FLUX PUREMA®				HIGH FLUX PUREMA®	MediSulfone®	
Membrane area (m²)	0,7	0,9	1,2	1,5	0,3	0,075	0,25
Fiber wall thickness (µm)	30				30	50	50
Housing inner diameter (mm)	38				24	14	24
Priming volume (ml)	49	57	75	87	21	5	19
Max transmembrane pressure	600 mmHg				600 mmHg		
Potting compound	Polyuteran				Polyuteran		
Housing and ferrule material	Polycarbonate (PC)				Copolyester		
Blood connectors	Twist Lock				Twist Lock	Male Luer lock	Twist Lock
Effluent connectors	Female Luer Lock				Female Luer Lock		
Ultrafiltration coefficient* (ml/h/mmHg)	21	27	36	45	9	1,3	5,5
Sieving coefficient	Albumin: about 0,002**				Albumin: about 0,002**	Albumin: 0,004***	
Cut-off	40.000 Daltons				40.000 Daltons	50.000 Daltons	
Sterilization	ETO - Beta rays				ETO - Beta rays		
Sterilization validity	3 years				3 years		
Product code ETO sterilized UDI-DI code	M03913 8033772326186	M03893 8033772326025	M03895 8033772326049	M03897 8033772326063	M03912 8033772326179	M03698 8033772322379	M03069 8033772320054
Product code Beta sterilized UDI-DI code	M03916 8033772326209	M03901 8033772326100	M03903 8033772326124	M03905 8033772326148	M03915 8033772326193	M03745 8033772323161	M03137 8033772320061
Quantity/box	6				10	15	10

Specifications are subject to change without notice. This document has no contractual value.

HANSEN CONNECTION

	DP60HE	DP120HE	DP150HE	DP190HE	DP230HE
Intended use	Adult				
Hollow fiber material	HIGH FLUX PUREMA®				
Membrane area (m²)	0,6	1,2	1,5	1,9	2,3
Fiber wall thickness (µm)	30				
Housing inner diameter (mm)	32		36	40	44
Priming volume (ml)	39	69	87	115	133
Max transmembrane pressure (mmHg)	600				
Potting compound	Polyurethane				
Housing and ferrule material	Polycarbonate (PC)				
Blood connectors	Twist Lock				
Effluent connectors	Hansen				
Ultrafiltration coefficient* (ml/h/mmHg)	18	36	45	57	69
Sieving coefficient	Albumin: about 0,002**				
Cut-off	40.000 Daltons				
Sterilization	ETO - Beta rays				
Sterilization validity	3 years				
Product code ETO sterilized UDI-DI code	M03892 8033772326018	M03894 8033772326032	M03896 8033772326056	M03898 8033772326070	M03899 8033772326087
Product code Beta sterilized UDI-DI code	M03900 8033772326094	M03902 8033772326117	M03904 8033772326131	M03906 8033772326155	M03907 8033772326162
Quantity/box	10	28		20	
DEHP free Latex free					

BRO-HMF-ENG_REV00 - CDP240725A

*Tests performed according to ISO8637 on bovine blood, Hct 32 ± 3%, TP = 60 ± 5 gr/L

** (QF = 30 ml/(min x m²))

*** Ref. ISO 8637 5.6.2

PUREMA® is a registered trademark of Membrana GmbH

MediSulfone® is a registered trademark of Medica SpA

DISTRUBITOR

MEDICA
SEAMLESS INNOVATION FOR LIFE QUALITY

MANUFACTURER
MEDICA S.p.A.
Via degli Artigiani, 7
41036, Medolla (MO), Italy

Phone +39 0535 51159
Fax + 39 0535 52605
info@medica-spa.com
www.medica-spa.com

CE 0476
CE 0476
UNI EN ISO 9001:2015
UNI CEI EN ISO 13485:2021